



### Sequences of Enhancer Elements (continued)

Vector	Sequence of Enhancer Element
pHTS-AP1	( <u>TGACTAA</u> ) <sub>6</sub>
pHTS-GAS	( <u>AGTTTCATATTACTCTAAATC</u> ) <sub>4</sub>
pHTS-ISRE	( <u>TAG TTT CAC TTT CCC</u> ) <sub>5</sub>
pHTS-CRE	(AGCC <u>TGACGTCA</u> GAG) <sub>4</sub>
pHTS-NFAT	( <u>GGAGGAAAACTGTTTCATACAGAAGGCGT</u> ) <sub>4</sub>
pHTS-NFκB	(T <u>GGGGACTTTC</u> GC) <sub>5</sub>

(Consensus sequences of the respective enhancer elements are underlined; subscripts indicate the number of tandem repeats of the sequence within the brackets in the vector.)

### Inserting new enhancer elements into pHTS-MCS

Cloning of your own *cis*-response elements: Complementary oligos of chosen *cis*-response elements may be synthesized and annealed with appropriate overhangs to be cloned in between Sma I and Sall of pHTS-MCS. For example, the following annealed dsDNA can be inserted into the Sall site of pHTS-MCS to create a reporter vector that is responsive to TGF-β.

5' -TCGATCTCAATCCACAATCTCGGAGTATGTCTAGACTGACAATG -3'  
 3' - AGAGTTAGGTGTTAGAGCCTCATACAGATCTGACTGTTACAGCT -5'

### Transfection and selection of mammalian cells

The provided plasmids are highly purified and ready for transfection into mammalian cells with any standard transfection methods such as LipofectAmine®, electroporation or Calcium Phosphate Precipitation. Refer to the manufacturers' instruction manuals for details of these transfection methods. Clones that harbor integrated pHTS vectors can be selected from the resulting transfected cells with hygromycin to establish stable cell lines.

### Luciferase Activity Assay

There're several commercial sources for reagents prepared for luciferase extraction and activity assay (e.g. Promega, Stratagene). The following protocol is provided for quick reference only.

1. Remove media from cell and rinse twice with PBS and remove residual PBS.
2. Add 1x Lysis Buffer (e.g. 400µl per well of a six-well plate, see below for buffer components). Incubate the plate for 15 minutes at room temperature (RT) with gentle rocking.
3. The lysates could be used in luciferase assay or be transferred to microcentrifuge tubes and stored at -80°C.
4. Mix 5-20µl of cell lysate with 100µl of 1X Luciferase Assay Buffer in an appropriate tube (e.g. Falcon® 2054 polystyrene tube). All reagents should be brought to RT for assay.
5. Measure the light emission with a luminometer with an integration time of 10-30 seconds.

Buffers for Luciferase Activity Assay: (Final concentrations are shown)

Lysis Buffer (5X)		Assay Buffer (1X)	
40mM	Tricine (pH7.8)	40mM	Tricine (pH7.8)
50mM	NaCl	0.5mM	ATP
2mM	EDTA	10mM	MgSO <sub>4</sub>
1mM	MgSO <sub>4</sub>	0.5mM	EDTA
5mM	DTT	10mM	DTT
1%	Triton® X-100	0.5mM	Coenzyme A
		0.5mM	luciferin

### Key References

1. Boulikas T. (1995) Phosphorylation of transcription factors and control of the cell cycle. *Critical Review in Eukaryotic Gene Expression* 5: 1-77
2. Decker, T. et al., 1991. Cytoplasmic activation of GAF, an IFN-γ-regulated DNA binding factor. *EMBO J.*, 10: 927-932
3. Fiering, S. et al., 1990. Single cell assay of a transcription factor reveals a threshold in transcription activated by signals emanating from the T-cell antigen receptor. *Genes Dev.*, 4: 1823-1834.
4. Galien R, Emanoil-Ravier R, Mercier G. (1994) Differential effects of c-jun and CREB on c-AMP response element activation by Ha-ras. *Oncogene* 9: 1101-1108
5. Karin, M. and Hunter, T. (1995) Transcriptional control by protein phosphorylation: signal transmission from the cell surface to the nucleus. *Curr. Biol.* 5: 747-757
6. Levy, DE., Kessler, DS, Pine, R, Reich, N. and Darnell, JE, Jr., 1998, Interferon-induced nuclear factors that bind a shared promoter element correlate with positive and negative transcriptional control. *Genes Devel.* 2: 383-393).
7. Treisman, R. (1996) Regulation of transcription by MAP kinase cascades. *Current Opinion in Cell Biol.* 8: 205-215
8. Wingender, E. (1990) Transcription regulating proteins and their recognition sequences. *Critical Review in Eukaryotic Gene Expression* 1: 11-48

### Notes

Falcon®, Geneticine®, LipofectAmine® and Triton® are registered trademarks of their respective owners.