

P53 ProtéGene™ Set

Catalog# P1035
 Lot# 171021

Materials Provided

1. pMEV-P53-WT (P1035a): 20 µg in 40 µl TE (pH7.5), 0.5 mg/ml.
2. pMEV-P53-DN (P1035b): 20 µg in 40 µl TE (pH7.5), 0.5 mg/ml.
3. Product Information Sheets.

Note: Individual plasmids can be ordered separately. Some plasmids are shipped as lyophilized pellet.

Receiving and Storage:

If received in liquid form, spin the vials briefly in a microcentrifuge to collect the contents. If received in lyophilized form, add 40 µl sterile DI water to the vial, mix thoroughly by vortex and then collect the contents by centrifuging the vials briefly in a microcentrifuge. Store the products at 2-8°C if used immediately and store at -20°C for extended storage.

Expression Vector:

pMEV-2HA (a): Cat# P1001a.

Affinity Tag:

N-terminal 2 x HA, a 9-aa peptide derived from influenza virus (MGYPYDVPDYAYPYDVPDYAGS...).

Prokaryotic Selection:

The kanamycin-resistance gene (aminoglycoside 3' phosphotransferase) expression cassette in the plasmids confers Kanamycin resistance to bacteria cells. Bacterial cells transformed with the plasmids should be maintained and grown in media containing 25-50µg/ml Kanamycin (e.g. #LK-1100, Prepared LB Agar plates, Biomyx, San Diego, California).

Eukaryotic Selection:

The neomycin resistance gene, driven by SV40 early promoter, confers G418 resistance to eukaryotic cells. Stable mammalian cell lines can be selected with G418.

Description of P53 and Mutants

The p53 protein is a nuclear phosphoprotein of 53 kDa. It is a tumor suppressor protein critical for regulation of the cell cycle in response to genotoxic insults such as radiation or chemicals resulting in DNA damages. The highly conserved vertebrate p53 gene is one of the most frequently mutated genes found in human cancers, being either lost or mutated in over half of all human tumors See Cox and Lane, 1995; Levine, 1997 for review). It has therefore become the center of intensive study ever since the link of p53 mutations with various human cancers was realized. The expression level of several important genes including p21, MDM2, GADD45 (CHOP), Bax and cyclin G was modulated by p53. Among many other roles, it is involved in arresting the cell cycle at G1 phase upon DNA damages (Cox and Lane, 1995; Levine, 1997).

The native p53 protein is a tetramer in solution and functions as a transcription factor that binds to specific DNA sequences. Many P53 mutations, especially those lie within the DNA-binding domain (so-called hot-spot mutations), show dominance over wild type protein. Replacement of valine 143, one of the hot-spot residues, results in a mutant (V143A) with a transdominant-negative effect on several activities of wild type p53. V143A is also a temperature-sensitive mutant: it shows DNA binding and transcriptional activity at 32°C, but it is inactive at 37°C (see Wong et al., 1999; Bartussek et al., 2002, for example).

Molecular Features of the Inserts:

Gene: *Homo sapiens* tumor protein p53 (Li-Fraumeni syndrome) (TP53)

GenBank Reference Sequence: NM_000546

5'-Cloning Site: Bam HI

5'-Junction Sequence(upper strand):

5'...tac gct gga tcc **ATG GAG GAG**...3'

3'-Cloning Site: Eco RI

3'-Junction Sequence (lower strand):

5'...acgcgtgaattc **TCA GTC TGA GTC** ...3'

P53 Protein Sequence

(393 amino acid residues. V143 is in bold and underlined.)

MEEPQSDPSVEPPLSQETFSDLWKLLENVLSPLPSQAMDDLMLSPDDI
 EQWFTEDPGDEAPRMEAPRVAAPAAPTAAAPAPAPSWPLSSSVPSQ
 KTYQGSYGRFLRGLHSGTAKSVTCTYSPALNKMFCQLAKTCTVQLWVDST
 PPPGTRVRAMAIYKQSQHMTEVVRRCPPHHERCSDSDGLAPPQHLIRVEGN
 LRVEYLDLDRNTFRHSVVVPEPEVGVSDCTTIHYNMNCSSCMGGMNRRP
 ILTIITLEDSSGNLLGRNSFEVRCACPRDRRTEENLRKKKGEPHHEL
 PGSTKRALPNNTSSSPQPKKPLDGEYFTLQIRGRERFEMFRELNEALEL
 KDAQAGKEPGGGRAHSSHLKSKKQGSTSRHKKLMFKTEGPDSD

P53 Nucleotide Sequence

(1182 bps. Codon for V143 is in bold and underlined)

```

1 ATGGAGGAGC CGCAGTCAGA TCCTAGCGTC GAGCCCCCTC TGAGTCAGGA
51 AACATTTTCA GACCTATGGA AACTACTTCC TGA AAAACAAC GTTCTGTCCC
101 CCTTGCCGTC CCAAGCAATG GATGATTTGA TGCTGTCCCC GGACGATATT
151 GAACAATGGT TCACTGAAGA CCCAGGTCCA GATGAAGTCC CCAGAATGCC
201 AGAGGCTGCT CCCCGCGTGG CCCCTGCACC AGCGACTCTT ACACCGCGCG
251 CCCCTGCACC AGCCCCCTCC TGGCCCCCTGT CATCTTCTGT CCCTTCCCAG
301 AAAACCTACC AGGGCAGCTA CGGTTCCCGT CTGGGCTTCT TGCATCTCTG
351 GACAGCCAAG TCTGTGACTT GCACGTACTC CCCTGCCTCT AACAAAGATG
401 TTTGCCAACT GGCCAAGACC TGCCCTGTG AGCTGTGGGT TGATTCACAC
451 CCCCGCCCGG GCACCCGCGT CCGCGCATG GCCATFACA AGCAGTCACA
501 GCACATGACG GAGGTTGTGA GGCCTGCACC CCACCATGAC CGCTGCTCAG
551 ATAGCGATGG TCTGGCCCTC CCTCAGCATC TTATCCGAGT GGAAGGAAAT
601 TTGCGTGTGG AGTATTTTGA TGACAGAAAC ACTTTTTCGAC ATAGTGTGGT
651 GTGCGCCTAT GAGCCGCTTG AGGTTGGCTC TGACTGTACC ACCATCCACT
701 ACAACTACAT GTGTAAACAGT TCCTGCATGG GCGGCATGAA CCGGAGGCCC
751 ATCCTCACCA TCATCACACT GGAAGACTCC AGTGGTAATC TACTGGGACG
801 GAACAGCTTT GAGGTGCGTG TTTGTGCTGT TCCTGGGAGA GACCCGCGCA
851 CAGAGGAAGA GAATCTCCCG AAGAAGGGGG AGCCTCACCA CGAGTCTGCC
901 CAGGGGAGCA CTAAGCGAGC ACTGCCCAAC AACACCAGCT CCTCTCCCA
951 GCCAAAGAAG AAACCACTGG ATGGAGAATA TTTCAACCTT CAGATCCCGT
1001 GCGGTGAGCG CTTGAGATG TCCGAGAGC TGAATGAGGC CTGGAAGATC
1051 AAGGATGCCC AGGCTGGGAA GGAGCCAGGG GGGAGCAGGG CTCACCTCAG
1101 CCACCTGAAG TCCAAAAGG GTCAGTCTAC CTCCGCCCAT AAAAAACTCA
1151 TGTTCAAGAC AGAAGGGCCT GACTCAGACT GA
  
```

Mutations:

pMEV-P53-WT (P1035a): No mutation

pMEV-P53-DN (P1035b): V143A (GTG to GCG)

Selected References:

- Tang C-H and Grimm E, 2004. Depletion of Endogenous Nitric Oxide Enhances Cisplatin-induced Apoptosis in a p53-dependent Manner in Melanoma Cell Lines. *J. Biol. Chem.* 279:288-298
- Bartussek C, Naumann U, Weller M. 2002. Tumour p53 mutations exhibit promoter selective dominance over wild type p53. *Oncogene.* 21:1641-8
- Cox, L.S. and Lane, D.P., 1995. Tumor suppressors, kinases and clamps: how p53 regulates the cell cycle in response to DNA damage. *BioEssays.* 17: 501-508
- Kern, S.E., Pietenpol, J.A., Thiagalingam, S., Seymour, A., Kinzler, K.W., and Vogelstein, B., 1992. Oncogenic forms of p53 inhibit p53-regulated gene expression. *Science.* 256: 827-830
- Levine, A., 1997. p53, the cellular gate keeper for growth and division. *Cell.* 88:323-331
- Wong KB, DeDecker BS, Freund SM, Proctor MR, Bycroft M, Fersht AR., 1999. Hot-spot mutants of p53 core domain evince characteristic local structural changes. *Proc Natl Acad Sci USA.* 96: 8438-42
- Zakut-Houri, R., Bienz-Tadmor, B, Givol, D. and Oren, M., 1985. Human p53 cellular tumor antigen: Cdna sequence and expression in COS cells. *The EMBO J.* 4: 1251-1255

Zambetti, G. P., Bargonetti, J., Walker, K., Prives, C. and Levine, A. J.,
1992, Wild-type p53 mediates positive regulation of gene expression
through a specific DNA sequence element. *Gene & Dev.* , 6: 1143-
1152